

INTRA - / EXTRAORAL

Intraorale Elastics - Latex

- Hergestellt aus chirurgischem Latex
- Zip-Patienten-Tüte mit Farb Code zur leichten Größenidentifizierung
- Standardpackung, 1.000 Elastics




Intraoral Elastics – Latex

- Made by surgical latex
- Force in ounces is measured at 3x passive diameter
- Simple to dispense and easy for patients to use
- Packaged in color-coded tamper-proof patient bags

VERPACKUNGSEINHEIT



100 Elastics pro Tüte



Standardpackung = 10 Tüten
100x10=1000 Elastics



Großpackung = 50 Tüten
100x50=5000 Elastics

Innendurchmesser Inner diameter	Leicht 2,5 oz (ROT) Light 2,5 oz (RED)	Mittel 3,5 oz (GRÜN) Medium 3,5 oz (GREEN)	Stark 4,5 oz (BLAU) Strong 4,5 oz (BLUE)	Extra stark 6 oz (BRAUN) Extra strong 6 oz (BROWN)
1/8" = 3,2 mm (A)	G731	G733	G735	G736
3/16" = 4,8 mm (B)	G741	G743	G745	G746
1/4" = 6,4 mm (C)	G751	G753	G755	G756
5/16" = 7,9 mm (D)	G761	G763	G765	G766
3/8" = 9,5 mm (E)	G771	G773	G775	G776

100 Elastics pro Tüte
1 Packung = 10 Tüten

100 elastics/bag
1 package = 10 bags = 1.000 Elastics

Intraorale Elastics - Latex

- Großpackung, 5.000 Elastics

Intraoral Elastics – Latex

- Jumbo package, 5.000 Elastics

Innendurchmesser Inner diameter	Leicht 2,5 oz (ROT) Light 2,5 oz (RED)	Mittel 3,5 oz (GRÜN) Medium 3,5 oz (GREEN)	Stark 4,5 oz (BLAU) Strong 4,5 oz (BLUE)	Extra stark 6 oz (BRAUN) Extra strong 6 oz (BROWN)
1/8" = 3,2 mm (A)	G731H	G733H	G735H	G736H
3/16" = 4,8 mm (B)	G741H	G743H	G745H	G746H
1/4" = 6,4 mm (C)	G751H	G753H	G755H	G756H
5/16" = 7,9 mm (D)	G761H	G763H	G765H	G766H
3/8" = 9,5 mm (E)	G771H	G773H	G775H	G776H

100 Elastics pro Tüte
1 Packung = 50 Tüten

100 elastics/bag
1 package = 50 bags = 5.000 Elastics