


BRACKETS



Brackets

Keramik Bracket Ultraclear

MBT

▲ OK / Upper			.018"	.018"	.022"	.022"
	Torque	Angulation	Left	Right	Left	Right
U1	17°	4°	AC-31-311-12	AC-31-311-11	AC32-2021	AC32-2011
U2	10°	8°	AC-31-311-21	AC-31-311-12	AC32-2022	AC32-2012
U3 w/hook	-7°	8°	AC-31-311-23HK	AC-31-311-13HK	AC32-2023H	AC32-2013H
U4	-7°	0°	AC-31-311-24	AC-31-311-14	AC32-2024	AC32-2014
U4 w/hook	-7°	0°	AC-31-311-24HK	AC-31-311-14HK	AC32-2024H	AC32-2014H
U5	-7°	0°	AC-31-311-25	AC-31-311-15	AC32-2025	AC32-2014
U5 w/hook	-7°	0°	AC-31-311-25HK	AC-31-311-15HK	AC32-2025H	AC32-2014H
▼ UK / Lower			.018"	.018"	.022"	.022"
	Torque	Angulation	Left	Right	Left	Right
U1	-6°	0°	AC-31-311-31	AC-31-311-41	AC32-2031	AC32-2041
U2	-6°	0°	AC-31-311-32	AC-31-311-42	AC32-2032	AC32-2042
U3 w/hook	-6°	3°	AC-31-311-33HK	AC-31-311-43HK	AC32-2033H	AC32-2043H
U4	-12°	2°	AC-31-311-34	AC-31-311-44	AC32-2034	AC32-2044
U4 w/hook	-12°	2°	AC-31-311-34HK	AC-31-311-44HK	AC32-2034H	AC32-2044H
U5	-17°	2°	AC-31-311-35	AC-31-311-45	AC32-2035	AC32-2044
U5 w/hook	-17°	2°	AC-31-311-35HK	AC-31-311-45HK	AC32-2035H	AC32-2044H

1 Packung = 10 Stück

1 Package = 10 pieces